



A549 Host Cell Proteins

Immunoenzymetric Assay for the Measurement of A549 Host Cell Proteins Catalog # F230

Intended Use

This kit is intended for use in determining the presence of protein contamination from the human cell line A549 in products manufactured by recombinant expression in A549 host cells. The kit is for Research and Manufacturing Use Only, and is not intended for diagnostic use in humans or animals.

Summary and Explanation

Recombinant expression by the A549 cell line can be used for cost effective production of a number of complex products. When applied with viral vectors, A549 cells have also shown utility as production systems for various gene therapy agents. Many of these recombinant products are intended for use as therapeutic agents in humans and animals and as such must be highly purified. The manufacturing and purification process of these products leaves the potential for contamination by host cell proteins (HCPs) from the A549 cell line. Such contamination can reduce the efficacy of the therapeutic agent and result in adverse toxic or immunological reactions and thus it is desirable to reduce HCP contamination to the lowest levels practical. Immunological methods using antibodies to HCPs such as Western Blot and ELISA are widely accepted due to their specificity and sensitivity. While Western blot is a powerful method aiding in the identity of HCPs, it suffers from a number of limitations. Western blot is a complex and technique dependent procedure requiring a subjective interpretation of results. Furthermore, it is essentially a qualitative method and does not lend itself to obtaining quantitative answers. The sensitivity of Western blot is severely limited by the volume of sample that can be tested and by interference from the presence of high concentrations of the intended product. While Western Blot may be able to detect HCPs in samples from upstream in the purification process it often lacks adequate sensitivity to detect HCPs in purified downstream and final product. The microtiter plate immunoenzymetric assay (ELISA) method employed in this kit overcomes the limitations of Western blots providing on the order of 100 fold better sensitivity. This simple to use, highly sensitive, objective, and semi-quantitative ELISA is a powerful method to aid in optimal purification process development, process control, and in routine quality control and product release testing. This kit is "generic" in the sense that it is intended to react with essentially all of the HCPs that could contaminate the product independent of the purification process. The antibodies have been generated against and affinity purified using a mild lysate washed of A549 cells to obtain HCPs typically encountered in your initial product recovery step. Special procedures were utilized in the

generation of these antibodies to insure that low molecular weight and less immunogenic contaminants as well as high molecular weight components would be represented. As such this kit can be used as a process development tool to monitor the optimal removal of host cell contaminants as well as in routine final product release testing. Because of the high sensitivity and broad reactivity of the antibodies, this generic kit has been successfully validated for testing of final product HCPs in various products regardless of growth and purification process. When the kit can be satisfactorily validated for your samples, the application of a more process specific assay is probably not necessary in that such an assay may only provide information redundant to this generic assay. However, if your validation studies indicate the antibodies in this kit are not sufficiently reactive with your process specific HCPs it may be desirable to also develop a more process specific ELISA. This later generation assay may require the use of a more specific and defined antisera. Alternatively, if the polyclonal antibody used in this kit provides sufficient sensitivity and broad antigenic reactivity, it may be possible to substitute the standards used in this kit for ones made from the contaminants that typically co-purify through your purification process and thus achieve better accuracy for process specific HCPs. The suitability of this kit for a given sample type and product must be determined and validated experimentally by each laboratory. The use of a process specific assay with more defined antigens and antibodies in theory may yield better sensitivity however such an assay runs the risk of being too specific in that it may fail to detect new or atypical contaminants that might result from some process irregularity or change. For this reason it is recommended that a broadly reactive "generic" host cell protein assay be used as part of the final product purity analysis even when a process specific assay is available.

Principle of the Procedure

The A549 HCP assay is a two-site immunoenzymetric assay. Samples containing A549 HCPs are reacted in microtiter strips coated with an affinity purified capture antibody. A second horseradish peroxidase (HRP) enzyme labeled anti-A549 antibody is reacted simultaneously resulting in the formation of a sandwich complex of solid phase antibody - HCP - enzyme labeled antibody. The microtiter strips are washed to remove any unbound reactants. The substrate tetramethyl benzidine (TMB) is then reacted. The amount of hydrolyzed substrate is read on a microtiter plate reader and is directly proportional to the concentration of A549 HCPs present.

Reagents & Materials Provided

Component	Product #
Anti-A549:HRP Affinity purified goat antibody conjugated to HRP in a protein matrix with preservative, 1x12mL	F231
Anti-A549 coated microtiter strips 12x8 well strips in a bag with desiccant	F232
A549 HCP Standards Solubilized A549 HCPs in a bovine serum albumin matrix with preservative. Standards at 0, 2, 8, 25, 75, & 200ng/mL, 1 mL/vial	F233
TMB Substrate 3,3',5,5' Tetramethyl benzidine, 1x12mL	F005
Wash Concentrate (20X) Tris buffered saline with preservative, 1x50mL	F004
Stop Solution 0.5N sulfuric acid, 1x12mL	F006

Storage & Stability

- * All reagents should be stored at 2°C to 8°C for stability until the expiration date printed.
- * The substrate reagent should not be used if its stopped absorbance at 450nm is greater than 0.1.
- * Reconstituted wash solution is stable until the expiration date of the kit.

Materials & Equipment Required But Not Provided

Microtiter plate reader spectrophotometer with dual wavelength capability at 450 & 650nm.
(If your plate reader does not provide dual wavelength analysis you may read at just the 450nm wavelength.)
Pipettors - 50µL and 100µL
Repeating or multichannel pipettor - 100µL
Microtiter plate rotator (150 - 200 rpm)
Sample Diluent (recommended Cat # F233A)
Distilled water
1 liter wash bottle for diluted wash solution

Precautions

For research or manufacturing use only.

Stop reagent is 0.5N H₂SO₄. Avoid contact with eyes, skin, and clothing. At the concentrations used in this kit none of the other reagents are believed to be harmful.

This kit should only be used by qualified technicians.

Preparation of Reagents

- * Bring all reagents to room temperature.
- * Dilute wash concentrate to 1 liter in distilled water, and label with kit lot and expiration date and store at 4°C.

Procedural Notes

1. Complete washing of the plates to remove excess unreacted reagents is essential to good assay reproducibility and sensitivity. We advise against the use of automated or other manual operated vacuum aspiration devices for washing plates as these may result in lower specific absorbances, higher non-specific absorbance, and more variable precision. The manual wash procedure described below generally provides lower backgrounds, higher specific absorbance, and better precision. If duplicate CVs are poor or if the absorbance of the "0" standard minus a substrate blank is greater than 0.250, evaluate plate washing procedure for proper performance.

2. High Dose Hook Effect may be observed in samples with very high concentrations of HCP. High Dose Hook Effect is due to insufficient excess of antibody for very high concentrations of HCPs present in the sample. Samples greater than 200 µg/mL may give absorbances less than the 200 ng/mL standard. It is also possible for samples to have certain HCPs in concentrations exceeding the amount of antibody for that particular HCP. In such cases the absorbance of the undiluted sample may be lower than the highest standard in the kit however these samples will fail to show acceptable dilutional recovery/parallelism as evidenced by an apparent increase in HCP concentration with increasing dilution. High Dose Hook is most likely to be encountered from samples early in the purification process. If a hook effect is possible, samples should also be assayed diluted. If the HCP concentration of the undiluted sample is less than the diluted sample this may be indicative of the hook effect. Such samples should be diluted at least to the minimum required dilution (MRD) where the dilution adjusted value remains essentially constant. The HCP value to be reported for such samples is the dilution corrected value at or greater than the established MRD. The diluent used should be compatible with accurate recovery. The preferred diluent is our Cat# F233A available in 100mL, 500mL, or 1 liter bottles. This is the same material used to prepare the kit standards. As an alternative, sample diluent Catalog # I028 should provide comparable recovery. Other prospective diluents must be tested for recovery by using them to dilute the 200ng/mL standard, as described in the "Limitations" section below.

Limitations

- * Before relying exclusively on this assay to detect host cell proteins, each laboratory should validate that the kit antibodies and assay procedure yield acceptable specificity, accuracy, and precision. A suggested protocol for this validation can be obtained by contacting our Technical Services Department or at our web site.
- * The standards used in this assay are comprised of A549 HCPs solubilized by mechanical disruption and detergent. Western blot analysis of the antibodies used in this kit demonstrates that they recognize the majority of distinct PAGE separated bands seen using a sensitive protein

staining method like silver stain or colloidal gold. Because the vast majority of HCPs will be conserved among all strains of A549 this kit should be adequately reactive to HCPs from your strain. Other clients have successfully validated this kit for their individual strains demonstrating acceptable specificity, accuracy, and sensitivity for process intermediate samples as well as final product. However, there can be no guarantee that this assay will detect all proteins or protein fragments from your process. In recognition of this potential limitation, we suggest reporting unknown samples in arbitrary dose units such as “ng/mL or parts per million of total immuno-reactive HCP equivalents”.

* Certain sample matrices may interfere in this assay. The standards used in this kit attempt to simulate typical sample protein and matrices. However the potential exists that the product protein or other components in the sample matrix may result in either positive or negative interference in this assay. It is advised to test all sample matrices for interference by diluting the 200ng/mL standard, 1 part to 4 parts of the matrix containing no or very low HCP contaminants. This diluted standard when assayed as an unknown should give a value of 30 to 50 ng/mL. Consult *Cygnus Technologies* Technical Service Department for advice on how to quantitate the assay in problematic matrices.

* Avoid the assay of samples containing sodium azide (NaN_3) which will destroy the HRP activity of the conjugate and could result in the under-estimation of HCP levels.

Assay Protocol

*The assay is very robust such that assay variables like incubation times, sample size, and other sequential incubation schemes can be altered to manipulate assay performance for more sensitivity, increased upper analytical range, or reduced sample matrix interference. Before modifying the protocol from what is recommended, users are advised to contact our technical services for input on the best way to achieve your desired goals.

* The protocol specifies the use of an approved microtiter plate shaker or rotator for the immunological step. If you do not have such a device it is possible to incubate the plate without shaking however it will be necessary to extend the first immunological incubation step by one hour for a total of 4 hours to achieve comparable results to the 3 hour shaking protocol. **Do not shake during the 30 minute substrate incubation step as this may result in higher backgrounds and worse precision.**

* Thorough washing is essential to proper performance of this assay. Automated plate washing systems or other vacuum aspiration devices are not recommended. The manual method described in the assay protocol is preferred for best precision, sensitivity and accuracy. A more detailed discussion of this procedure can be obtained from our Technical Services Department or on our web site.

* Bring all reagents to room temperature. Set-up plate spectrophotometer to read dual wavelength at 450nm for the test wavelength and 650nm for the reference.

* Thorough washing is essential to proper performance of this assay.

* All standards, controls, and samples should be assayed in duplicate.

* Maintain a repetitive timing sequence from well to well for all assay steps to insure that all incubation times are the same for each well.

* Make a work list for each assay to identify the location of each standard, control, and sample.

* If the substrate has a distinct blue color prior to the assay it may have been contaminated. If the absorbance of 100 μL of substrate plus 100 μL of stop against a water blank is greater than 0.1 it may be necessary to obtain new substrate or the sensitivity of the assay may be compromised.

* Strips should be read within 30 minutes after adding stop solution since color will fade over time.

Assay Protocol

1. Pipette 50 μL of standards, controls and samples into wells indicated on work list.

2. Pipette 100 μL of anti-A549:HRP (#F231) into each well.

3. Cover & incubate on rotator at ~ 180rpm for 3 hours at room temperature, 24°C \pm 4°.

4. Aspirate, or dump the contents to waste. Blot and vigorously bang out residual liquid over absorbent paper. Wash generously with diluted wash solution by flooding the wells with solution from a squirt bottle or by pipeting in ~350 μL . Dump and bang again. Repeat for a total of 4 washes. Wipe off any liquid from the bottom outside of the microtiter wells as any residue can interfere in the reading step.

5. Pipette 100 μL of substrate (#F005).

6. Incubate for 30 minutes. Do Not Shake!

7. Pipette 100 μL of Stop Solution (#F006).

8. Read absorbance at 450/650nm blanking on the Zero standard.

Calculation of Results

The standards may be used to construct a standard curve with values reported in ng/mL “total immuno-reactive HCP equivalents” (See Limitations section above). This data reduction may be performed through computer methods using curve fitting routines such as point to point, cubic spline, or 4 parameter logistic fit. **Do not use linear regression analysis to interpolate values for samples as this may lead to significant inaccuracies!** Data may also be manually reduced by plotting the absorbance values of the standard on the y-axis versus concentration on the x-axis and drawing a smooth point to point line. Absorbances of samples are then interpolated from this standard curve.

Example Data

Well #	Contents	Abs. at 450nm	Mean Abs.	ng/mL HCP equivs.
1A	Zero Std	0.000		
1B	Zero Std	0.002	0.001	
1C	2ng/mL	0.020		
1D	2ng/mL	0.030	0.025	
1E	8ng/mL	0.094		
1F	8ng/mL	0.084	0.088	
1G	25ng/mL	0.240		
1H	25ng/mL	0.246	0.243	
2A	75ng/mL	0.688		
2B	75ng/mL	0.708	0.698	
2C	200ng/ml	1.908		
2D	200ng/ml	1.892	1.900	
2E	sample A	1.990		
2F	sample A	1.960	1.975	>200ng
2G	sample B	0.085		
2H	sample B	0.089	0.087	8.0ng

Quality Control

- Precision on duplicate samples should yield average % coefficients of variation of less than 10% for samples greater than 8 ng/mL and < 200 ng/mL. CVs for samples < 8ng/mL may be greater than 10%.

- For optimal performance the absorbance of the substrate when blanked against water should be < 0.1.

- It is recommended that each laboratory assay appropriate quality control samples in each run to insure that all reagents and procedures are correct.

Performance Characteristics

Cygnus Technologies has validated this assay by conventional criteria as indicated below. This validation is generic in nature and is intended to supplement, but not

replace certain user and product specific validation, which should be performed by each laboratory.

Sensitivity

The lower limit of detection (**LOD**) defined as that concentration corresponding to a signal two standard deviations above the mean of the zero standard is 0.5ng/mL.

The lower limit of quantitation (**LOQ**) defined as that concentration where concentration coefficients of variation (CVs) are <20% is 1.6ng/mL.

Precision

Both intra (n=20 replicates) and inter-assay (n=5 assays) precision were determined on 3 pools with low (~9ng/mL), medium (~25ng/mL), and high concentrations (~75ng/mL). The % CV is the standard deviation divided by the mean multiplied by 100.

Pool	Intra assay CV	Inter assay CV
Low	6.6	7.4
Medium	2.9	4.4
High	4.9	6.5

Recovery/ Interference Studies

A number of buffer matrices were evaluated by adding known amounts of the A549 HCP preparation used to make the standards in this kit. Because this assay is designed to minimize matrix interference most of these buffers yielded acceptable recovery defined as between 80-120%. In general extremes in pH (<5.0 and >8.5) as well as certain detergents can cause under-recovery. High concentrations of your product protein can also inhibit binding of the antibodies to A549 HCPs and thus result in an underestimation of true HCP levels. When such product or sample matrix interference can be demonstrated, it may be necessary to dilute or otherwise process your samples to achieve acceptable recovery. Recovery experiments can be simply performed by diluting some of the 200ng/mL standard provided with this kit into the sample type in question and assaying this ‘spiked’ sample for the expected HCP level. Consult *Cygnus Technologies* Technical Services Department for advice on how to conduct and interpret recovery experiments.

Specificity/Cross-Reactivity

Because of the broad reactivity of the antibodies used in this kit, the possibility exists for some immunological cross-reactivity to related or conserved proteins in your samples. Such true immunological cross-reactivity can manifest itself as either a false increase in apparent HCP levels when little or none exists, or a false decrease in true HCP levels if the cross-reactant is in sufficient excess to compete with HCP for a limited amount of antibody. It is also a possibility that some components in your samples, while not immunologically cross-reactive, may give a non-specific increase in assay signal that could be interpreted as either A549 HCPs or immunological cross-reactivity. Where practical you should evaluate components in your

samples for positive interferences in the form of either cross reactivity or increase non-specific binding. Contact *Cygnus Technologies* Technical Service Department for advice on how to conduct experiments to determine cross-reactivity and non-specific, positive interferences.

Hook Capacity

Increasing concentrations of HCPs > 200 ng/mL were assayed as unknowns. The hook capacity, defined as that concentration which will give an absorbance reading less than the 200 ng/mL standard was ~200 µg/mL.

Ordering Information/ Customer Service

To place an order or to obtain additional product information contact *Cygnus Technologies* Customer Support:

Tel: 910-454-9442 Fax: 910-454-9443
Email: cygnustec@aol.com
Web: www.cygnustechnologies.com
4701 Southport Supply Road SE, Suite 7
Southport, NC 28461 USA