



Cygnus Technologies, Inc.

4705 Southport Supply Rd. SE, Suite 208
Southport, NC 28461 USA
TEL: 910-454-9442 FAX: 910-454-9443
E-mail: cygnustec@aol.com

Material Safety Data Sheet

Catalog Number: IM40

Product Name: Human Immunoglobulin Class Screening Kit

The below information is believed to be correct but may not be all inclusive. It is intended to be used as a guide. You, as the user, should independently verify the suitability and safety of the product under your conditions of use. Cygnus Technologies shall not be liable or responsible for the use of this information or the material supplied and thus the user shall assume all risk in the use of the product. Disposal of hazardous material may be subject to local, state or federal regulations.

Description of Contents:

Goat, anti-hIgA:Horseradish Peroxidase, Product # IM140

Affinity purified goat antibody conjugated to horseradish peroxidase
ProClin 300 preservative
NaCl, (sodium chloride)
Tris {hydroxymethyl}-aminomethane

Goat, anti-hIgD:Horseradish Peroxidase, Product # IM141

Affinity purified goat antibody conjugated to horseradish peroxidase
ProClin 300 preservative
NaCl, (sodium chloride)
Tris {hydroxymethyl}-aminomethane

Goat, anti-hIgE:Horseradish Peroxidase, Product # IM142

Affinity purified goat antibody conjugated to horseradish peroxidase
ProClin 300 preservative
NaCl, (sodium chloride)
Tris {hydroxymethyl}-aminomethane

Goat, anti-hIgG, gamma chain specific:Horseradish Peroxidase, Product # IM143

Affinity purified goat antibody conjugated to horseradish peroxidase
ProClin 300 preservative
NaCl, (sodium chloride)
Tris {hydroxymethyl}-aminomethane

Goat, anti-hIgM, mu chain specific:Horseradish Peroxidase, Product # IM144

Affinity purified goat antibody conjugated to horseradish peroxidase
ProClin 300 preservative
NaCl, (sodium chloride)
Tris {hydroxymethyl}-aminomethane

Conjugate Diluent, Product # IM103

Tris {hydroxymethyl}-aminomethane
ProClin 300 preservative
Fish Gelatin

Hazardous Ingredients:

Only the ingredients listed below have any known toxicities. For the other materials shown above no biological hazard has been established. At the concentrations used the hazards are minimal, however this kit should only be used by experienced laboratory personnel trained in the proper techniques and general laboratory safety. The following procedures apply for all components of this kit.

General Toxicological Information - Inhalation of any aerosols from these components may induce sensitization and cause allergic reactions in sensitized individuals.

First Aid Measures -

Skin Contact: Wash skin with soap and water. If symptoms occur consult a doctor.

Eye Contact: Rinse eyes with water for 10 minutes and contact a doctor.

Ingestion: Rinse mouth and throat with water and consult a doctor.

Inhalation: Exit from exposure area. If symptoms of irritation, sensitization or if breathing difficulties occur, contact a doctor.

Fire fighting measures: No special requirement. Water or foam are suitable. There are no special exposure hazards.

Accidental release procedures: Spilled material should be absorbed into paper or cloth towels with the area washed with water. Avoid splashing or aerosolization. Ensure sufficient ventilation.

Handling and storage: Store at 2°C to 8°C. Avoid aerosolization.

Exposure controls and personal protection: Wear impermeable gloves, protective eye wear.

Hazardous Chemicals:

1. Tris{hydroxymethyl}-aminomethane

- Irritant to eyes, respiratory system, and skin

2. NaCl (sodium chloride)

- May be toxic if ingested in large quantities

- May cause irritation when contacted with eyes or skin

- Avoid inhalation

3. Proclin 300 Preservative

- Irritant to eyes, respiratory system, and skin